

# UFI FILTERS DISPOSABLE MEDICAL MASK FOR ADULTS


## THREE LAYERS OF PROTECTION

|                       |                         |
|-----------------------|-------------------------|
| FILTRATION EFFICIENCY | $\geq 95\%$             |
| DIFFERENTIAL PRESSURE | $< 40 \text{ Pa/cm}^2$  |
| MICROBIAL CLEANLINESS | $\leq 30 \text{ cfu/g}$ |



### 1 WATER-RESISTANT SPUNBOND

Filtration of particles of diameter  $3.0 \pm 0.3 \mu\text{m}$



### 3 HYGROSCOPIC SPUNBOND

Soft, comfortable and breathable

### 2 MELTBLOWN

Double strengthen, prevent missing dust



**UFI FILTERS**  
CHOSEN BY THE BEST

[www.ufifilters.com](http://www.ufifilters.com)

# UFI FILTERS DISPOSABLE MEDICAL MASK FOR ADULTS

## CHARACTERISTICS

- + Three-layers disposable medical mask with ear holder and adjustable nose clip.
- + Produced in accordance to EN 14683:2019 standard Type 1.
- + Distributed by UFI Filters Shanghai.
- + The masks are supplied in bag to preserve the cleanliness and contamination.
- + Must be stored in well-ventilated, dark and dry environment, away from fire and pollution.
- + Recommended hours of use: not more than 4 hours.
- + Shelf Life: 2 years.

NOSE CLIP FIXED BY  
ULTRASONIC WELDING



EAR HOLDER FIXED BY  
ULTRASONIC WELDING

## 20PCS – TRANSPARENT SINGLE PACKAGE

|  |            |   |           |
|--|------------|---|-----------|
|  | LIGHT BLUE |  | 35.011.71 |
|  | BLACK      |  | 35.011.72 |



2020 – Copyright © UFI FILTERS spa



**UFI FILTERS**  
**CHOSEN BY THE BEST**

[www.ufifilters.com](http://www.ufifilters.com)