

UFI FILTERS DISPOSABLE MEDICAL MASK FOR KIDS


THREE LAYERS OF PROTECTION

FILTRATION EFFICIENCY	≥ 95%
DIFFERENTIAL PRESSURE	< 40 Pa/cm ²
MICROBIAL CLEANLINESS	≤ 30 cfu/g



1 WATER-RESISTANT SPUNBOND

Filtration of particles of diameter $3.0 \pm 0.3 \mu\text{m}$



3 HYGROSCOPIC SPUNBOND

Soft, comfortable and breathable

2 MELTBLOWN

Double strengthened, prevent missing dust



UFI FILTERS
CHOSEN BY THE BEST

www.ufifilters.com

UFI FILTERS DISPOSABLE MEDICAL MASK FOR KIDS

CHARACTERISTICS

- + Three-layers disposable medical mask with ear holder and adjustable nose clip.
- + Produced in accordance to EN 14683:2019 standard Type 1.
- + Distributed by UFI Filters Shanghai.
- + The masks are supplied in bag to preserve the cleanliness and contamination.
- + Must be stored in well-ventilated, dark and dry environment, away from fire and pollution.
- + Recommended hours of use: not more than 4 hours.
- + Shelf Life: 2 years.

NOSE CLIP FIXED BY
ULTRASONIC WELDING



EAR HOLDER FIXED BY
ULTRASONIC WELDING

20PCS – TRANSPARENT
SINGLE PACKAGE



 LIGHT
BLUE



35.011.70

2020 – Copyright © UFI FILTERS spa



UFI FILTERS
CHOSEN BY THE BEST

www.ufifilters.com