


# UFI FILTERS MASK KN95 MASK DEFENDER

## FIVE LAYERS OF PROTECTION

FILTRATION EFFICIENCY	≥ 95%	@0.3 µm NaCl, 85l/min
INHALATION RESISTANCE	≤ 350pa	@85l/min
EXHALATION RESISTANCE	≤ 250pa	@85l/min




**UFI FILTERS**  
CHOSEN BY THE BEST

[www.ufifilters.com](http://www.ufifilters.com)

# UFI FILTERS MASK KN95 MASK DEFENDER

## CHARACTERISTICS

- + Five-layers face mask with ear holder and adjustable nose clip. Ultrasonically welded.
- + Produced in accordance to GB 2626-2006 KN95.
- + The production line is in UFI Filters Shanghai.
- + Masks are supplied in single plastic bags to preserve the cleanliness and contamination - 20 pieces per box.
- + Must be stored in well-ventilated, dark and dry environment, away from fire and pollution.
- + Recommended storage conditions: temperature -20C°~ +30C° and relative humidity less than 85%.
- + Shelf Life: 3 years.



**20PCS – TRANSPARENT  
SINGLE PACKAGE**

<input type="checkbox"/>	<b>WHITE</b>		<b>35.004.70</b>
<input type="checkbox"/>	<b>BLACK</b>		<b>35.004.71</b>



2020 – Copyright © UFI FILTERS spa



**UFI FILTERS**  
**CHOSEN BY THE BEST**

[www.ufifilters.com](http://www.ufifilters.com)