


UFI FILTERS MASK FFP2 MASK DEFENDER

FOUR LAYERS OF PROTECTION

FILTRATION EFFICIENCY	≥ 98%	@ 0.6µm, 95l/min
INHALATION RESISTANCE	≤ 2.4mbar	@ 95l/min
EXHALATION RESISTANCE	≤ 3.0mbar	@ 160l/min




UFI FILTERS
CHOSEN BY THE BEST

www.ufifilters.com

UFI FILTERS MASK FFP2 MASK DEFENDER

CHARACTERISTICS

- + Four-layers face mask with ear holder and adjustable nose clip. Ultrasonically welded.
- + Produced in accordance to EN 149:2001+ A1:2009 standard.
- + The production line is in UFI Filters Shanghai.
- + Masks are supplied in single plastic bags to preserve the cleanliness and contamination - 20 pieces per box.
- + Must be stored in well-ventilated, dark and dry environment, away from fire and pollution.
- + Recommended storage conditions: temperature -20C°~ +30C° and relative humidity less than 85%.
- + Shelf Life: 3 years.



20PCS – TRANSPARENT
SINGLE PACKAGE

☐ WHITE
☒ 35.005.70



2024 – Copyright © UFI FILTERS spa



UFI FILTERS
CHOSEN BY THE BEST

www.ufifilters.com